

Product Datasheet

Mouse TNFSF14 ELISA Kit (orb437549)

Description	Mouse Tumor Necrosis Factor Ligand Superfamily, Member 14 (TNFSF14) Sandwich ELISA Kit
Reactivity	Mouse
Range	15.625-1000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	Q9QYH9
Sensitivity	6.3pg/mL
Entrez	50930
Expiration Date	Please enquire.