

Product Datasheet

Human TNFSF9 ELISA Kit (orb437552)

Description	Human Tumor Necrosis Factor Ligand Superfamily, Member 9 (TNFSF9) Sandwich ELISA Kit
Reactivity	Human
Range	15.62-1000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P41273
Sensitivity	5.35pg/mL
Entrez	8744
Expiration Date	Please enquire.