

Product Datasheet

Equine VEGFA ELISA Kit (orb437569)

Description	Equine Vascular Endothelial Growth Factor A (VEGFA) Sandwich ELISA Kit
Reactivity	Equine
Range	31.25-2000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	Q9GKR0
Sensitivity	12.2pg/mL
Entrez	100033839
Expiration Date	Please enquire.