

Product Datasheet

Porcine VEGFA ELISA Kit (orb437572)

Description	Porcine Vascular Endothelial Growth Factor A (VEGFA) Sandwich ELISA Kit
Reactivity	Porcine
Range	15.625-1000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	P49151
Sensitivity	6.3pg/mL
Entrez	397157
Expiration Date	Please enquire.