

Product Datasheet

Human VEGFC ELISA Kit (orb437579)

Description	Human Vascular Endothelial Growth Factor C (VEGFC) Sandwich ELISA Kit
Reactivity	Human
Range	62.5-4000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, platelet-poor plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	serum, platelet-poor plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	P49767
Sensitivity	23.7pg/mL
Entrez	7424
Expiration Date	Please enquire.