

Product Datasheet

Bovine VEGFD ELISA Kit (orb437581)

Description	Bovine Vascular Endothelial Growth Factor D (VEGFD) Sandwich ELISA Kit
Reactivity	Bovine
Range	31.25-2000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	A7MBB6
Sensitivity	10.2pg/mL
Entrez	286799
Expiration Date	Please enquire.