

Product Datasheet

Human VEGFR1 ELISA Kit (orb437585)

Description	Human Vascular Endothelial Growth Factor Receptor 1 (VEGFR1) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, saliva, cell lysates, follicular fluid and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, saliva, cell lysates, follicular fluid and other biological fluids.
Uniprot ID	P17948
Sensitivity	0.066ng/mL
Entrez	2321
Expiration Date	Please enquire.