

Product Datasheet

Porcine VEGFR1 ELISA Kit (orb437587)

Description	Porcine Vascular Endothelial Growth Factor Receptor 1 (VEGFR1) Sandwich ELISA Kit
Reactivity	Porcine
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	F1RSU5
Sensitivity	0.115ng/mL
Entrez	396763
Expiration Date	Please enquire.