

Product Datasheet

CD262 antibody (PE) (orb43759)

Description

Mouse Monoclonal to CD262.

Reactivity

Human

Conjugation

PE

Tested

FC

Applications
Immunogen

Recombinant fusion protein of human IgG heavy chain and extracellular domain of DR5.

Target

CD262

Preservatives

Stabilizing phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide

Concentration

0.1 mg/ml

Storage

Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.

Note

For research use only

Application notes

Flow cytometry: Recommended dilution: 1-5 µg/ml; positive control: JURKAT human peripheral blood leukemia T cell line.

Isotype

Mouse IgG1

Clonality

Monoclonal

Purity

Purified antibody is conjugated with R-phycoerythrin (PE) under optimum conditions. Unconjugated antibody and free fluorochrome are removed by size-exclusion chromatography.

Uniprot ID

[O14763](#)

Entrez


[8795](#)

Dilution Range

Flow cytometry: Recommended dilution: 1-5 µg/ml; positive control: JURKAT human peripheral blood leukemia T cell line.

Expiration Date

12 months from date of receipt.



Flow
cytometry
analysis
(surface
stainin...