

Product Datasheet

Porcine VEGFR2 ELISA Kit (orb437591)

Description	Porcine Vascular Endothelial Growth Factor Receptor 2 (VEGFR2) Sandwich ELISA Kit
Reactivity	Porcine
Range	1.56-100ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	K7GNF0
Sensitivity	0.62ng/mL
Entrez	397311
Expiration Date	Please enquire.