

Product Datasheet

Rat VEGFR2 ELISA Kit (orb437592)

Description	Rat Vascular Endothelial Growth Factor Receptor 2 (VEGFR2) Sandwich ELISA Kit
Reactivity	Rat
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	008775
Sensitivity	0.063ng/mL
Entrez	25589
Expiration Date	Please enquire.