

Product Datasheet

Human VEGFR3 ELISA Kit (orb437593)

Description	Human Vascular Endothelial Growth Factor Receptor 3 (VEGFR3) Sandwich ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P35916
Sensitivity	0.125ng/mL
Entrez	2324
Expiration Date	Please enquire.