

Product Datasheet

Human VGF ELISA Kit (orb437595)

Description	Human VGF Nerve Growth Factor Inducible (VGF) Sandwich ELISA Kit
Reactivity	Human
Range	78.125-5000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	O15240
Sensitivity	32pg/mL
Entrez	7425
Expiration Date	Please enquire.