

Product Datasheet

Human IL12B ELISA Kit (orb437602)

Description	Human Interleukin 12B (IL12B) Sandwich ELISA Kit
Reactivity	Human
Range	31.25-2000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	P29460
Sensitivity	11.84pg/mL
Entrez	3593
Expiration Date	Please enquire.