

Product Datasheet

Human a1ACT ELISA Kit (orb437611)

Description	Human Alpha-1-Antichymotrypsin (a1ACT) Sandwich ELISA Kit
Reactivity	Human
Range	3.125-200ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P01011
Sensitivity	1.42ng/mL
Entrez	12
Expiration Date	Please enquire.