

Product Datasheet

Rat a2M ELISA Kit (orb437626)

Description	Rat Alpha-2-Macroglobulin (a2M) Sandwich ELISA Kit
Reactivity	Rat
Range	15.62-1000ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, urine and other biological fluids.
Sample Types	serum, plasma, urine and other biological fluids.
Uniprot ID	P06238
Sensitivity	8.7ng/mL
Entrez	100911545
Expiration Date	Please enquire.