

## Product Datasheet

### Human AASDHPPT ELISA Kit (orb437631)

|                          |                                                                                               |
|--------------------------|-----------------------------------------------------------------------------------------------|
| <b>Description</b>       | Human Aminoadipate Semialdehyde Phosphopantetheinyl Transferase (AASDHPPT) Sandwich ELISA Kit |
| <b>Reactivity</b>        | Human                                                                                         |
| <b>Range</b>             | 0.156-10ng/mL                                                                                 |
| <b>Form/Appearance</b>   | Traditional                                                                                   |
| <b>Storage</b>           | All reagents should be kept according to the instructions on the packaging upon receipt.      |
| <b>Note</b>              | For research use only                                                                         |
| <b>Application notes</b> | tissue homogenates and other biological fluids.                                               |
| <b>Sample Types</b>      | tissue homogenates and other biological fluids.                                               |
| <b>Uniprot ID</b>        | <b>Q9NRN7</b>                                                                                 |
| <b>Sensitivity</b>       | 0.054ng/mL                                                                                    |
| <b>Entrez</b>            | <b>60496</b>                                                                                  |
| <b>Expiration Date</b>   | Please enquire.                                                                               |