

Product Datasheet

Rat Ab1-40 ELISA Kit (orb437634)

Description	Rat Amyloid Beta Peptide 1-40 (Ab1-40) competitive inhibition ELISA Kit
Reactivity	Rat
Range	31.25-2000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Sensitivity	11.6pg/mL
Expiration Date	Please enquire.