

Product Datasheet

Human ABCC10 ELISA Kit (orb437642)

Description	Human ATP Binding Cassette Transporter C10 (ABCC10) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q5T3U5
Sensitivity	0.061ng/mL
Entrez	89845
Expiration Date	Please enquire.