

Product Datasheet

Human ABCC13 ELISA Kit (orb437645)

Description	Human ATP Binding Cassette Transporter C13 (ABCC13) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q9NSE7
Sensitivity	0.057ng/mL
Entrez	150000
Expiration Date	Please enquire.