

Product Datasheet

Human ABCC9 ELISA Kit (orb437653)

Description	Human ATP Binding Cassette Transporter C9 (ABCC9) Sandwich ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	O60706
Sensitivity	0.132ng/mL
Entrez	10060
Expiration Date	Please enquire.