

Product Datasheet

Mouse ABCG2 ELISA Kit (orb437656)

Description	Mouse ATP Binding Cassette Transporter G2 (ABCG2) Sandwich ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	Q7TMS5
Sensitivity	0.055ng/mL
Entrez	26357
Expiration Date	Please enquire.