

Product Datasheet

Human ACACa ELISA Kit (orb437658)

Description	Human Acetyl Coenzyme A Carboxylase Alpha (ACACa) Sandwich ELISA Kit
Reactivity	Human
Range	3.125-200ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q13085
Sensitivity	1.09ng/mL
Entrez	31
Expiration Date	Please enquire.