



Product Datasheet

CD45 antibody (orb43766)

Description	Rat monoclonal antibody to CD45
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	FC, ICC, IP, WB
Immunogen	Murine peripheral blood leukocytes
Target	CD45
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: Recommended dilution: 2 µg/ml; positive control: murine peripheral blood leukocytes.
Isotype	Rat IgG
Clonality	Monoclonal
Purity	Purified by protein-A affinity chromatography.
Uniprot ID	Q8C6Q7
Entrez	19264
Dilution Range	Flow cytometry: Recommended dilution: 2 µg/ml; positive control: murine peripheral blood leukocytes.
Expiration Date	12 months from date of receipt.



Separation of murine splenocytes stained...



Flow cytometry surface staining pattern ...