

Product Datasheet

Human ACE2 ELISA Kit (orb437663)

Description	Human Angiotensin I Converting Enzyme 2 (ACE2) Sandwich ELISA Kit
Reactivity	Human
Range	15.625-1000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q9BYF1
Sensitivity	7.7pg/mL
Entrez	59272
Expiration Date	Please enquire.