

Product Datasheet

Mouse ACE2 ELISA Kit (orb437664)

Description	Mouse Angiotensin I Converting Enzyme 2 (ACE2) Sandwich ELISA Kit
Reactivity	Mouse
Range	62.5-4000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	Q8R0I0
Sensitivity	24.4pg/mL
Entrez	70008
Expiration Date	Please enquire.