

Product Datasheet

Bovine ACE ELISA Kit (orb437665)

Description	Bovine Angiotensin I Converting Enzyme (ACE) Sandwich ELISA Kit
Reactivity	Bovine
Range	6.25-400ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q9BDG1
Sensitivity	2.15ng/mL
Entrez	509484
Expiration Date	Please enquire.