

Product Datasheet

Porcine ACE ELISA Kit (orb437668)

Description	Porcine Angiotensin I Converting Enzyme (ACE) Sandwich ELISA Kit
Reactivity	Porcine
Range	31.25-2000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	A1YQ85
Sensitivity	10.7pg/mL
Entrez	613133
Expiration Date	Please enquire.