

Product Datasheet

Rat ACE ELISA Kit (orb437669)

Description	Rat Angiotensin I Converting Enzyme (ACE) Sandwich ELISA Kit
Reactivity	Rat
Range	3.125-200ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, cerebrospinal fluid and other biological fluids.
Sample Types	serum, plasma, cerebrospinal fluid and other biological fluids.
Uniprot ID	P47820
Sensitivity	1.16ng/mL
Entrez	24310
Expiration Date	Please enquire.