


Product Datasheet

CD45 antibody (PerCP) (orb43769)

Description	Rat monoclonal antibody against CD45 conjugated to PerCP.
Reactivity	Mouse
Conjugation	PerCP
Tested Applications	FC
Immunogen	Murine peripheral blood leukocytes
Target	CD45
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	0.5 mg/ml
Storage	Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: Recommended dilution: 2-4 µg/ml; positive control: murine peripheral blood leukocytes.
Isotype	Rat IgG
Clonality	Monoclonal
Purity	Purified antibody is conjugated with activated Peridinin-Chlorophyll Protein (PerCP) under optimum conditions and unconjugated antibody and free fluorochrome are removed by size-exclusion chromatography.
Uniprot ID	Q8C6Q7
Entrez	19264
Dilution Range	Flow cytometry: Recommended dilution: 2-4 µg/ml; positive control: murine peripheral blood leukocytes.
Expiration Date	12 months from date of receipt.



Separation of murine splenocytes stained...