

Product Datasheet

Rat ACTH ELISA Kit (orb437693)

Description	Rat Adrenocorticotrophic Hormone (ACTH) competitive inhibition ELISA Kit
Reactivity	Rat
Range	15.625-1000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Sensitivity	4.76pg/mL
Expiration Date	Please enquire.