

Product Datasheet

Human ADAM10 ELISA Kit (orb437708)

Description	Human A Disintegrin And Metalloprotease 10 (ADAM10) Sandwich ELISA Kit
Reactivity	Human
Range	62.5-4000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	O14672
Sensitivity	24.6pg/mL
Entrez	102
Expiration Date	Please enquire.