

Product Datasheet

Rat ADAM10 ELISA Kit (orb437709)

Description	Rat A Disintegrin And Metalloprotease 10 (ADAM10) Sandwich ELISA Kit
Reactivity	Rat
Range	78.125-5000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	Q10743
Sensitivity	27pg/mL
Entrez	29650
Expiration Date	Please enquire.