

Product Datasheet

Human ADAM19 ELISA Kit (orb437714)

Description	Human A Disintegrin And Metalloprotease 19 (ADAM19) Sandwich ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	Q9H013
Sensitivity	0.108ng/mL
Entrez	8728
Expiration Date	Please enquire.