

Product Datasheet

Human ADAM28 ELISA Kit (orb437716)

Description	Human A Disintegrin And Metalloprotease 28 (ADAM28) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	Q9UKQ2
Sensitivity	0.061ng/mL
Entrez	10863
Expiration Date	Please enquire.