

Product Datasheet

Human ADAM8 ELISA Kit (orb437718)

Description	Human A Disintegrin And Metalloprotease 8 (ADAM8) Sandwich ELISA Kit
Reactivity	Human
Range	6.25-400pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	P78325
Sensitivity	2.35pg/mL
Entrez	101
Expiration Date	Please enquire.