

Product Datasheet

Human ADAMTS4 ELISA Kit (orb437724)

Description	Human A Disintegrin And Metalloproteinase With Thrombospondin 4 (ADAMTS4) Sandwich ELISA Kit
Reactivity	Human
Range	62.5-4000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	075173
Sensitivity	26.3pg/mL
Entrez	9507
Expiration Date	Please enquire.