

Product Datasheet

Rat ADAMTS7 ELISA Kit (orb437729)

Description	Rat A Disintegrin And Metalloproteinase With Thrombospondin 7 (ADAMTS7) Sandwich ELISA Kit
Reactivity	Rat
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q1EHB3
Sensitivity	0.058ng/mL
Entrez	315879
Expiration Date	Please enquire.