

Product Datasheet

Rat ADCY10 ELISA Kit (orb437732)

Description	Rat Adenylate Cyclase 10, Soluble (ADCY10) Sandwich ELISA Kit
Reactivity	Rat
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	Q9Z286
Sensitivity	0.125ng/mL
Entrez	59320
Expiration Date	Please enquire.