

Product Datasheet

Mouse ADCYAP1 ELISA Kit (orb437737)

Description	Mouse Adenylate Cyclase Activating Polypeptide 1, Pituitary (ADCYAP1) Sandwich ELISA Kit
Reactivity	Mouse
Range	78.125-5000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	070176
Sensitivity	28pg/mL
Entrez	11516
Expiration Date	Please enquire.