

Product Datasheet

Rat ADCYAP1 ELISA Kit (orb437738)

Description	Rat Adenylate Cyclase Activating Polypeptide 1, Pituitary (ADCYAP1) Sandwich ELISA Kit
Reactivity	Rat
Range	6.25-400pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P13589
Sensitivity	3.12pg/mL
Entrez	24166
Expiration Date	Please enquire.