

Product Datasheet

Rat ADH1 ELISA Kit (orb437741)

Description	Rat Alcohol Dehydrogenase 1 (ADH1) Sandwich ELISA Kit
Reactivity	Rat
Range	1.56-100ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P06757
Sensitivity	0.59ng/mL
Entrez	24172
Expiration Date	Please enquire.