

Product Datasheet

Human ADH2 ELISA Kit (orb437742)

Description	Human Alcohol Dehydrogenase 2 (ADH2) Sandwich ELISA Kit
Reactivity	Human
Range	15.625-1000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P00325
Sensitivity	6.7pg/mL
Entrez	125
Expiration Date	Please enquire.