

Product Datasheet

Bovine ADM ELISA Kit (orb437750)

Description	Bovine Adrenomedullin (ADM) Sandwich ELISA Kit
Reactivity	Bovine
Range	15.625-1000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	O62827
Sensitivity	6.4pg/mL
Entrez	280713
Expiration Date	Please enquire.