

Product Datasheet

Rat ADM ELISA Kit (orb437755)

Description	Rat Adrenomedullin (ADM) Sandwich ELISA Kit
Reactivity	Rat
Range	15.625-1000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P43145
Sensitivity	3.87pg/mL
Entrez	25026
Expiration Date	Please enquire.