

Product Datasheet

Human ADORA2b ELISA Kit (orb437759)

Description	Human Adenosine A2b Receptor (ADORA2b) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	P29275
Sensitivity	0.061ng/mL
Entrez	136
Expiration Date	Please enquire.