

Product Datasheet

Mouse ADORA2b ELISA Kit (orb437760)

Description	Mouse Adenosine A2b Receptor (ADORA2b) Sandwich ELISA Kit
Reactivity	Mouse
Range	78.125-5000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q60614
Sensitivity	30pg/mL
Entrez	11541
Expiration Date	Please enquire.