

Product Datasheet

General ADP ELISA Kit (orb437761)

Description	General Adenosine Diphosphate (ADP) competitive inhibition ELISA Kit
Reactivity	All
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Sensitivity	0.101ng/mL
Expiration Date	Please enquire.