

Product Datasheet

Human ADRa2C ELISA Kit (orb437763)

Description	Human Adrenergic Receptor Alpha 2C (ADRa2C) Sandwich ELISA Kit
Reactivity	Human
Range	0.625-40ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	P18825
Sensitivity	0.237ng/mL
Entrez	152
Expiration Date	Please enquire.