

Product Datasheet

Rat ADRb1 ELISA Kit (orb437766)

Description	Rat Adrenergic Receptor Beta 1 (ADRb1) Sandwich ELISA Kit
Reactivity	Rat
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	P18090
Sensitivity	0.055ng/mL
Entrez	24925
Expiration Date	Please enquire.